



# Perspectives on Synchrotron Radiation Research and the AOFSRR

### Richard Garrett Senior Advisor, Synchrotron Science, ANSTO



2013 Cheiron School

# Asia/Oceania Forum for Synchrotron Radiation Research

#### **AOFSRR Objective**

The objective of the AOFSRR is to <u>encourage regional collaboration</u> in, and to <u>promote the advancement of</u>, <u>synchrotron radiation research</u> and related subjects in Asia and Oceania.

#### **Specific Activities:**

- (1) The annual workshop and Cheiron School, and organization of other scientific collaboration meetings;
- (2) Exchange of information of facilities and user groups;
- (3) Provision of a framework for cooperative activities;
- (4) Any activities that promote and expand the role of synchrotron light source facilities and synchrotron based research in the Asia – Oceania region.



## **European Union**

Large: ESRF (+PETRA) Medium: Diamond, Soleil, SLS Soft X 3rg gen: Elettra, Max, BESSY II ... Next generation: FLASH, European XFEL, FERMI, PSI..



**ESRF** 

Diamond

**FLASH** 

# The Americas

APS

CLS

LCLS (Stanford)

Large: APS Medium: CLS, NSLS II, SSRL Soft X 3rd Gen: ALS 2<sup>nd</sup> gen: NSLS, CHESS, Brazil, Alladin... Next generation: LCLS, JLab, Cornell ERL(?)



#### Photon Factory

## Asia Oceania & the AOFSRR





Siam



- Facilities equal or better than Europe & USA
- Many bi-lateral agreements between facilities
- Few relationships between user communities
- No real regional organisation







PF,

ore ...

### **Associate Members:**

- Malaysia
- New Zealand
- Vietnam

### **Members:**

- Australia
- China
- India
- Japan
- South Korea
- Singapore
- Taiwan
- Thailand



# Asia/Oceania Forum for Synchrotron Radiation Research



synchrotron Radiation

#### The 1st AOFSRR Summer School Cherron School SPring-8, Japan

September 10th – 20th 2007 Organizer: AOFSRR, RIKEN, Spring-8, JASJRI, KEK-P

The Ordern School's main an available study of the design of the second state of the second s

er Natt Helpel Stadiogen American Provident of ALVSIN, Japon Helpel Stadio Stadiogen Der Provident of ALVSIN, Japon Helpel Maschi Takata (INEXINJASRI String 4, Apopt

Const Multier Richter Grenz (Anstalla) Ourras Stimsfram (Installa) In See Kaffkranz) Bibber Maser (Stiggawic) Right Nachar (Colub Managang Jatinot) (Takal) Bingli Anti (Jana) Bertain Kernel (Anstalla) Stell, J. Chang (Januar) Maethor Reco (Kernel)

end Departing Menther Mondel Takan (EUX 2014) Thing 41 Manaye Sanati (2014) ENAng 40 Kanali Sentrasch (2014) ENAng 40 Hannet Thianna (1/ASRI Strings 40 Hannet (2014) ENAng 40 Hannet Salama (1/ASRI Shrings 40 Vanhanet Salama (1/ASRI Shrings 40 Salari Salama (1/ASRI Shrings 40 Salari Salama (1/ASRI Shrings 40

ternalistal Advisory Beard "To be source on Enzyold Sternates (RECEVERS) Vers Kern (ASERDagen) A.G. "Hill" Stoling (ESBET-Ennet) Soling (March (Darmond VK) Verheim Konet (MS Darm)

http://cheiron2007.spring8.or.jp

### AOFSRR 2012 Bright light for better life

6<sup>th</sup> Asia-Oceania Forum for Synchrotron Radiation Research And 4<sup>th</sup> SLRI Annual User Meeting

> August 8-12, 2012 Imperial Queen's Park, Bangkok, Thailand

# **AOFSRR Activities**



# Annual Workshop

Year	Host		
2006	Tsukuba, Japan		
2007	Hsinchu, Taiwan		
2008	Melbourne, Australia		
2009	Shanghai, China		
2010	Pohang, South Korea		
2011/12	Bangkok, Thailand		
2013	Himeji, Japan		
2014	Hsinchu, Taiwan		

### Cheiron School: Always SPring-8 !!



# **User Community Networking**

- The AOF annual workshop
- Cheiron School
- Open access to facilities
- Special Access
- Multi-nation scientific collaborations
- Regional accelerator school
- Other workshops



# Promote Synchrotron Research

- Nations can build new communities at other facilities
  - Australian soft X-ray program at NSRRC and
  - NSRRC hard X-ray program at SPring-8
- Promote SR research in non-member nations in the A-O region
- Assist SR science in developing nations
  - Cheiron School
  - Assistance to attend conferences
  - Visiting scientists hosted at SR facilities
  - Work with other organisations (IUCr etc)



# In the Future it is <u>your</u> AOF: How Should it Develop ?



## Synchrotron Radiation



X25 wiggler beam, NSLS



## Outline

- What is a synchrotron?
- How is the light produced & what are its characteristics?
- Brief Basics of Synchrotron Beamlines
- Some Applications
- The Future (is here already): "Next Generation Sources"
- Some Cool Stuff



### **Development of SR Sources**



Courtesy M. Ree, PLS







#### Section of the Australian Synchrotron







S [m]

Generation of Synchrotron Radiation: Radiation from Accelerating Charge



Low energy electrons: Radiation in all directions Example: Radio waves from a transmitter.



High energy (relativistic) electrons:

Radiation pattern swept into a narrow cone in the forward direction = High brightness!



#### **Third Generation Sources: Undulator Insertion Devices**









AST 210/EECS 213 Univ. California, Berkeley

## **Beamline Design Goals**

- Deliver the required X-ray beam to the experiment:
  - Energy and bandwidth
  - Spot size
  - Divergence/convergence
- Preserve source characteristics eg intensity, coherence
- Optimise signal / background
- Be very stable and reproducible, in position, intensity and energy
- Be safe to operate
- Be user friendly to operate
- Achieve all the above within a reasonable budget !





#### Hard X-ray Beamline: Si crystal monochromator E > 4 keV



#### Soft X-ray Beamline: Grating monochromator E < 2 keV



## **Mirrors for Synchrotron Beamlines**

- Deflection
- Focusing
- Harmonic Rejection
- Power Reduction





## Critical Angle/Reflectivity with Energy: Rhodium Coated Mirror Example



Harder X-rays need more grazing angles and longer mirrors:

2 mm high beam needs:

≤ 10 cm mirror at 1 keV

≥ 80 cm mirror at 20 keV



## An example beamline: the AS Xray Absorption Spectroscopy Beamline



## Characteristics of Synchrotron Radiation

- ✓ High brightness/flux
- ✓ Wide energy spectrum
- ✓ Plane polarised
- ✓ Pulsed





## Unique Characteristics of Synchrotron Radiation

- Extremely high brightness. Modern synchrotron sources are about 10 billion times as intense as a laboratory X-ray generator: dilute samples; fast measurements; trace elements;
  - Low divergence: high intensity can be focussed onto tiny samples: Microscopies
- •Wide X-ray energy spectrum:
  - the optimum X-ray energy to be chosen for each experiment;
  - > X-ray spectroscopies are possible eg EXAFS
- Polarisation: various dichroisms; magnetic imaging; molecular orientation;
- Time structure: time of flight and very fast timing.



### Why Brightness is Important



**Diamond Anvil Cell** ~ 100 μm sample chamber

Brilliance of a synchrotron (flux/source size/divergence) gives more usable x-rays or effective intensity



## X-rays and their Interaction with Matter



X-ray Fluorescence → trace element analysis

Transmitted Photons:
→ Imaging
Absorption Spectroscopy
→ Chemical information



#### X-ray Diffraction -> Structure

	Synchrotron	Proton	Electron Microscope	SIMS	Neutron
Sensitivity	~	•	×	~	×
Sub micron	$\checkmark$	•	$\checkmark$	×	×
Chemical Information	~	×	•	•	×
In-situ	~	•	×	×	<b>~</b>
Atomic Structure	<b>√</b>	×	~	×	<b>~</b>



## Sometimes High Intensity = Better Data Synchrotron Powder XRD



Multi-phase ceramic:  $\alpha$  Al<sub>2</sub>O<sub>3</sub>, ZrO<sub>2</sub>, MgO-Al<sub>2</sub>O<sub>3</sub> (spinel). Top synchrotron data; Bottom: lab data.

B. H. Oconnor, A. van Riessen, J. Carter, G. Burton, R. F. Garrett and D. J. Cookson, J. American Chemical Soc. 80 (1997) 1373



### X-ray Diffraction at a Synchrotron



## Most Significant Macromolecular Structures Solved using Synchrotron X-rays







The fusion protein from Newcastle Disease Virus. ©Structure M. Lawrence, CSIRO



Virus



Neuraminidase: Colman & Varghese, CSIRO

# **C**<sub>12</sub> а Fo b $\mathbf{F}_{\mathbf{1}}$ δ α

ATP Synthase: a Molecular Motor

H. Wang and G. Oster (1998). Nature 396:279-282.

John Walker won the 1997 Nobel Chemistry prize for solving the F1 catalytic domain using synchrotron radiation at Daresbury, UK.



#### Atomic structure informs biological function



## Broad Energy Spectrum: SR Only Spectroscopies eg Xray Absorption Spectroscopy

#### XANES: near edge structure

Sensitive to chemical environment of absorbing element.

Often different valence states have markedly different XANES spectra.



benign) and Cr VI, a known carcinogen.

EXAFS: extended structure to ~1 keV above an absorption edge

Nearest neighbour atomic distances, coordination etc. Crystals not required: disordered systems like solution species can be measured.



Amorphous GaAs EXAFS and Fourier transform.



### Near Edge Spectroscopy: Chemical Sensitivity



Carbon K-edge Spectra





Pt spectrum located in a tumour cell Hambley, U Syd





Courtesy H Ade, NC State Univ.

### **Micro-imaging and XANES from Mineral Sands**

- Aim: to image and analyse chemical state of radioactive trace elements Zircon mineral grains
- X-ray Absorption Near Edge Spectroscopy (XANES) determined uranium is present as U<sub>IV</sub>
- Chemical state is needed to design process to remove these elements



Zr K-edge 18.0 keV U L3-edge 17.2 keV



Concentration ranges:

U 13 ppm to 33 ppm Th 3 ppm to 11 ppm Pb 3 to 4 ppm





Radiography













Imaging

## Some Imaging Needs Focusing optics

Reflective (Kirkpatrick-Baez mirrors)typical ~1 μmHigh efficiency, achromatic, limited to ~10 nm

Diffractive (Fresnel zone plates) Moderate efficiency, limited to ~10 nm

Refractive (compound refractive lenses) 10s μm - ~ 50 nm Low efficiency, highly chromatic, aberrations



typical ~ 100 nm

FOCUS





## Contrast mechanisms in x-ray imaging

- Absorption measure electron density; can be element specific
- Fluorescence measure elemental distribution
- Spectroscopy extract chemical state, spin state
- Diffraction reveal structure, strain, magnetism, charge...
- Phase measure real part of refractive index

In general with X-rays:

- Natural sample contrast is possible; staining not required
- Image structure of thick samples, sectioning not required
- More penetrating, less damage, less charging than with electrons



# Phase Contrast

Refractive index:

for X-rays it is less than 1 by about 1 part in a million

$$n = 1 - \delta - i\beta = 1 - \frac{r_e}{2\pi}\lambda^2 \sum_i n_i f_i(0)$$

- Absorption contrast: sensitive to *Im(n)*
- Phase contrast: sensitive to Re(n)
- At high X-ray energies, phase contrast wins





## Phase Contrast

2. Diffragionio En Plance Contraistg (DEI)



AS Medical Beamline Wiggler Source Phase contrast simulations (Can recover phase shift)









# **Diffraction Enhanced Imaging**

- Edge/density gradient sensitive
- Move on rocking curve to change contrast





## X-ray phase imaging: Biology and Materials





Synchrotron



Rob Lewis, Monash University/ Energy = 20keV

Synchrotron: propagation phase contrast



DEI "Dark field" phase image of bonded aluminium sheets @ 33 keV Dots are bubbles in the epoxy bond.

Stevenson, Garrett, Hyodo et.al.



# **Two Famous Microscopes**



## First: ALS XM-1 Microscope



## Nanoscale 3-D biotomography

Mother daughter yeast cells just before separation



2-D slice from 3-D Tomogram. Images every 2°, 180° data set, several minutes.  $\Delta r = 45$  nm

Color coding identifies subcellular components by their x-ray absorption coefficients



Courtesy of Carolyn Larabell, UCSF/LBNL

### Second: The X-ray Fluorescence Microscope Beamline at the AS



Conventional Fluorescence Microscope: APS 2ID-D

## Advanced fluorescence detector at the AS

### **Annular geometry**

Maximises solid angle, sample @ 90°
384 Si pixel detector array (BNL, Siddons etal)

• No constraint on lateral sample size and scan range



### + Parallel data processing

- CSIRO: HYMOD2 pipelined, parallel processor (Ryan etal)
- Whole XRF spectrum acquired and analysed in real time

### + Fast Scanning Stage

Data acquired "on the fly"
milli-second dwell times
cf 1 second or greater normally

= New micro-XRF capability at the AS X-ray Fluorescence Microscope beamline



#### XFM image definition (number of pixels) limited by dwell time

## Long dwell → Low Image Definition

- ~1 s / pixel (for readout of 1-16 detector spectra)
- 1.3 hours  $\rightarrow$  67 x 67 pixels





#### Short dwell → Good image definition

- 32 ms / pixel
- 1.3 hours → 375 x 375 pixels (30 fold increase)
- New Maia-32 prototype detector, NSLS X27A



Courtesy C. Ryan, CSIRO

## New Sources:

**XFELs** 



### Linac-Based Free Electron Laser Self-Amplified Spontaneous Emission (SASE)







<u>Remarkable Features of XFEL producing  $\lambda < 0.1$  nm X-Rays</u>

- O High Peak Brilliance
- <u>Narrow Pulse Width</u>
- **O** High Degree of Coherence



#### SACLA 1<sup>st</sup> beamline: 90m Undulator



### Coherent Diffractive Imaging: no lens, no crystal



R. Gerchberg and W. Saxton, Optik 35, 237 (1972); J.R. Fienup, Appl. Opt. 21, 2758 (1982)

## Single Shot Imaging at the FLASH Soft X-ray FEL





A coherent diffraction pattern of the object recorded from a single 25-femtosecond FEL pulse.



Reconstructed image: no signs of damage caused by the pulse.

#### H. Chapman et al., Nature Physics 2, 839 (2006)

### **XFEL Holy Grail?: Single Molecule Imaging**









Neutze etal, Nature 406, (2000) 752-757

### Nano Crystallography - seems most promising so far

pН

3.0

3.5

4.0

3.0

4.5

8.0

6.0

[protein] (mg/ml)

3d protein nano-crystals

avoid damage ⇒ pulse < 20 fs





2D membrane protein crystals





R. Abela, PSI

# **Two Cool Examples to Finish**



## **Exploring Cultural Heritage**

Portrait of a Woman By Edgar Degas 1876-1880

Courtesy David Thurrowgood (National Gallery of Victoria)





## **Exploring Our Cultural Heritage**





 $3.7 \text{ ms per } 60 \times 60 \ \mu\text{m2} \ \text{pixel}$ 

32 megapixel dataset



## Phase Contrast CT: Fossils in Amber





Courtesy P Tafforeau, ESRF